



Product Datasheet

XKR8 antibody (orb45542)


Description	XKR8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human XK-related protein 8 protein (69-158AA)
Target	XKR8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H6D3
Dilution Range	WB:1:500-5000, IHC-P:1:1000-2000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human li...