

Product Datasheet

BCL11B antibody (Biotin) (orb45569)

Description	BCL11B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human B-cell lymphoma/leukemia 11B protein (500-550AA)
Target	BCL11B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9C0K0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.