



Product Datasheet

FLRT3 antibody (orb45586)

Description	FLRT3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Leucine-rich repeat transmembrane protein FLRT3 protein (209-264AA)
Target	FLRT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZU0
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart(lan...



Immunofluorescence analysis of HepG2 cel...