


Product Datasheet

AFF4 antibody (orb45650)

Description	AFF4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human AF4/FMR2 family member 4 protein (77-334AA)
Target	AFF4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UHB7
Dilution Range	WB:1:2000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver tis...