

Product Datasheet

PPM1H antibody (Biotin) (orb45673)

Description	PPM1H antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Protein phosphatase 1H protein (73-289AA)
Target	PPM1H
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ULR3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.