

Product Datasheet

LRCH1 antibody (Biotin) (orb45705)

Description	LRCH1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Leucine-rich repeat and calponin homology domain-containing protein 1 protein (1-105AA)
Target	LRCH1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2L9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.