

Product Datasheet

Eukaryotic aspartyl protease antibody (Biotin)
[Out of stock] (orb457432)

Description

Rabbit polyclonal antibody

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide
derived from Oryza nivara
Eukaryotic aspartyl protease

Preservatives

0.01M TBS(pH7.4) with 1%
rAlbumin, 0.03% Proclin300 and
50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year.
Avoid repeated freeze/thaw
cycles. The lyophilized antibody is
stable at room temperature for at
least one month and for greater
than a year when kept at -20°C.
When reconstituted in sterile pH
7.4 0.01M PBS or diluent of
antibody the antibody is stable
for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

52kDa

Dilution Range

ELISA: 1:100-1000, IHC-P: 1:50-
200, IF: 1:50-200

Expiration Date

12 months from date of receipt.