

Product Datasheet

GAD65 antibody (Biotin) (orb457445)

Description	GAD65 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human GAD65
Target	GAD65
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q05329
Dilution Range	WB=1:500-2000, ICC=1:100, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.