



Product Datasheet

C5AR1 antibody (orb45762)

Description	C5AR1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human C5a anaphylatoxin chemotactic receptor 1 protein (1-36AA)
Target	C5AR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21730
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...



Western blot analysis of PC-3 whole cel...