



Product Datasheet

CD40 antibody (orb45774)

Description	CD40 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 5 protein (21-193AA)
Target	CD40
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25942
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse lung(lane...



Immunofluorescence analysis of HeLa cell...