

Product Datasheet

HBsAg antibody (Biotin) (orb457802)

Description	HBsAg antibody (Biotin)
Species/Host	Goat
Reactivity	Human, Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant HBsAg protein
Target	HBsAg
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q89714
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.