

Product Datasheet

EFEMP1 antibody (Biotin) (orb45817)

Description	EFEMP1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human EGF-containing fibulin-like extracellular matrix protein 1 protein (27-184AA)
Target	EFEMP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12805
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.