



Product Datasheet

EIF4A3 antibody (orb45822)


Description	EIF4A3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Eukaryotic initiation factor 4A-III protein (6-166AA)
Target	EIF4A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38919
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunofluorescence analysis of MCF-7 cells.



Western blot analysis of A549 whole cell lysate.