

Product Datasheet

NT 3 antibody (Biotin) (orb459565)

Description	NT 3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Neurotrophin 3
Target	Neurotrophin 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20783
Dilution Range	IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.