


Product Datasheet

MPI antibody (orb45986)

Description	MPI antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Mannose-6-phosphate isomerase protein (213-350AA)
Target	MPI
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P34949
Dilution Range	IF:1:20-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...