

Product Datasheet

Keratocan antibody (Biotin) (orb460102)

Description	Keratocan antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Keratocan
Target	Keratocan
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60938
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.