

Product Datasheet

Transglutaminase 2 antibody (Biotin) (orb460116)

Description	Terminally differentiating mammalian epidermal cells acquire an insoluble, 10 to 20 nm thick protein...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Transglutaminase 2: 351-450aa/687aa
Target	Transglutaminase 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	77 kDa
Uniprot ID	P21980
Expiration Date	12 months from date of receipt.