

Product Datasheet

**SLC27A4/FATP4 antibody (Biotin) [Out of stock]
(orb460394)**

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Conjugation	Biotin
Immunogen	KLH conjugated synthetic peptide derived from human SLC27A4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72kDa
Uniprot ID	Q6P1M0, Q91VE0
Dilution Range	WB=1:50-200 ELISA=1:100-1000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.