

Product Datasheet

Beta galactosidase antibody (Biotin) [Out of stock] (orb460664)

Description	Beta galactosidase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human Beta galactosidase
Target	Beta galactosidase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71kDa
Uniprot ID	P16278
Dilution Range	WB: 1:500-1000, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.