

Product Datasheet

SH3BP5 antibody (Biotin) (orb46094)

Description	SH3BP5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human SH3 domain-binding protein 5 protein (49-455AA)
Target	SH3BP5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60239
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.