

Product Datasheet

TRX antibody (Biotin) (orb460946)

Description	TRX antibody (Biotin)
Species/Host	Mouse
Reactivity	All, Human, Porcine
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant human Thioredoxin
Target	Trx-Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	12kDa
Uniprot ID	P10599
Dilution Range	WB=1:1000-5000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.