

Product Datasheet

Tubulin beta 3 antibody (Biotin) (orb460949)

Description	Tubulin beta 3 antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TUBB3
Target	TUBB3 (Neuronal Marker)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q13509
Dilution Range	WB=1:500-5000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.