

Product Datasheet

Aquaporin 4 antibody (Biotin) (orb460988)

Description	Aquaporin 4 antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	Recombinant human Aquaporin 4 Protein
Target	Aquaporin 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P55087
Dilution Range	IHC-P=1:100-500,IHC-F=1:100- 500,ICC=1:100-500,ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.