

Product Datasheet

Dopamine D2 receptor antibody (Biotin) [Out of stock] (orb460995)

Description	Dopamine D2 receptor antibody (Biotin) [Out of stock]
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Dopamine D2 Receptor
Target	Dopamine D2 Receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14416
Dilution Range	WB: 1:500-2000, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.