

Product Datasheet

SMS antibody (Biotin) (orb46122)

Description	SMS antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Spermine synthase protein (16-191AA)
Target	SMS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P52788
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.