

Product Datasheet

Human CD45-FITC antibody (FITC) [Out of stock] (orb461283)

Description	Mouse monoclonal antibody against Human CD45-FITC conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Recombinant human CD45
Target	PTPRC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	143 kDa
Uniprot ID	P08575
Dilution Range	Flow-Cyt=20ul/Test
Expiration Date	12 months from date of receipt.