

Product Datasheet

HLA-DR-FITC antibody (FITC) (orb461293)

Description	Mouse monoclonal antibody against HLA-DR-FITC conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Recombinant human HLA-DR
Target	HLA-DRA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
Uniprot ID	P01903
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.