

## Product Datasheet

### Cytomegalovirus UL23 antibody (FITC) [Out of stock] (orb461472)

|                        |                                                                                                                                                                                                                                                                                                                               |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Rabbit polyclonal antibody against Cytomegalovirus UL23 conjugated to FITC                                                                                                                                                                                                                                                    |
| <b>Species/Host</b>    | Rabbit                                                                                                                                                                                                                                                                                                                        |
| <b>Conjugation</b>     | FITC                                                                                                                                                                                                                                                                                                                          |
| <b>Immunogen</b>       | KLH conjugated synthetic peptide derived from HCMV UL23 C-terminus                                                                                                                                                                                                                                                            |
| <b>Preservatives</b>   | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                                                                                                                                                                                                                                         |
| <b>Form/Appearance</b> | Lyophilized or Liquid                                                                                                                                                                                                                                                                                                         |
| <b>Concentration</b>   | - 100 µg (in 100 µl): 1 mg/ml                                                                                                                                                                                                                                                                                                 |
| <b>Storage</b>         | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                                                                                         |
| <b>Isotype</b>         | IgG                                                                                                                                                                                                                                                                                                                           |
| <b>Clonality</b>       | Polyclonal                                                                                                                                                                                                                                                                                                                    |
| <b>MW</b>              | 33kDa                                                                                                                                                                                                                                                                                                                         |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                                                                                                                                               |