

Product Datasheet

Beta-Actin (Loading Control) antibody (FITC) (orb461483)

Description	Beta-Actin (Loading Control) antibody (FITC)
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	Synthetic MAP peptide derived from human beta-Actin
Target	beta-Actin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60709
Dilution Range	Flow-Cyt=1µg/Test,ICC=1:100
Expiration Date	12 months from date of receipt.