

Product Datasheet

CHRM4/mAChR M4 antibody (FITC) (orb461669)

Description	CHRM4/mAChR M4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CHRM4/mAChR M4
Target	CHRM4
Preservatives	PBS, Glycerol, rAlbumin
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53kDa
Uniprot ID	P08173
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.