

Product Datasheet

**GRIK2 + GRIK3 antibody (FITC) [Out of stock]
(orb461682)**

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Conjugation	FITC
Immunogen	KLH conjugated synthetic peptide derived from human GRIK1 + GRIK2 + GRIK3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	99kDa
Uniprot ID	Q13002 , P42260 , P39087
Dilution Range	ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.