

Product Datasheet

CEA antibody (FITC) (orb461795)

Description	Mouse monoclonal antibody against CEA conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from Human CEA/CD66e/CEACAM5
Target	CEACAM5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	150-200kDa
Uniprot ID	P06731
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.