



Product Datasheet

TUBA4A antibody (orb46227)


Description	TUBA4A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Tubulin alpha-4A chain protein (20-183AA)
Target	TUBA4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68366
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3(lane 1)...



Immunohistochemical staining of human CO...



Immunofluorescence analysis of HepG2 cel...