

Product Datasheet

GSDML/Gasdermin like antibody (FITC) [Out of stock] (orb462292)

Description

Rabbit polyclonal antibody against

FITC

Species/Host

Rabbit

Conjugation

FITC

Immunogen

KLH conjugated synthetic peptide derived from human GSDML/Gasdermin like

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

47kDa

Uniprot ID
Q8TAX9
Dilution Range

ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.