


Product Datasheet

SLC2A2 antibody (orb46259)

Description	SLC2A2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Solute carrier family 2, facilitated glucose transporter member 2 protein (67-174AA)
Target	SLC2A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11168
Dilution Range	IF:1:20-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...