

Product Datasheet

MEIKIN antibody (Biotin) (orb46290)

Description	MEIKIN antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Meiosis-specific kinetochore protein (18-91AA)
Target	MEIKIN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0A087WXM9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.