

Product Datasheet

Gapdh1 antibody (orb46299)

Description	Gapdh1 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Drosophila melanogaster Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-332AA)
Target	Gapdh1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07486
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.