



## Product Datasheet

### PRKAA1 antibody (orb46339)


|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PRKAA1 antibody                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                            |
| <b>Reactivity</b>          | Human, Mouse                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                                                |
| <b>Immunogen</b>           | Recombinant Human 5\\'-AMP-activated protein kinase catalytic subunit alpha-1 protein (325-543AA) |
| <b>Target</b>              | PRKAA1                                                                                            |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | IgG                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q13131</a>                                                                            |
| <b>Dilution Range</b>      | WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |



Immunohistochemical staining of human ov...



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human lu...