


Product Datasheet

PLEKHA7 antibody (orb46395)

Description	PLEKHA7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Pleckstrin homology domain-containing family A member 7 protein (294-654AA)
Target	PLEKHA7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IQ23
Dilution Range	IF:1:20-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...