

Product Datasheet

Alpha Tubulin (Acetyl Lys40) antibody (FITC) [Out of stock] (orb464108)

Description

Mouse monoclonal antibody

Species/Host

Mouse

Reactivity

Human, Mouse, Zebrafish

Conjugation

FITC

Immunogen

KLH conjugated synthetic peptide derived from human alpha Tubulin (Acetyl Lys40)

Target

TUBA4A

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

MW

50kDa

Uniprot ID
P68366
Expiration Date

12 months from date of receipt.