

Product Datasheet

Bcl-2 antibody (FITC) (orb464158)

Description	Bcl-2 antibody (FITC)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	FITC
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human Bcl-2
Target	Bcl-2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P49950
Expiration Date	12 months from date of receipt.