

## Product Datasheet

### Dystrobrevin Alpha antibody (FITC) (orb464290)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Dystrobrevin Alpha antibody (FITC)                                               |
| <b>Species/Host</b>        | Rabbit                                                                           |
| <b>Reactivity</b>          | Human, Mouse                                                                     |
| <b>Conjugation</b>         | FITC                                                                             |
| <b>Tested Applications</b> | ICC, IF                                                                          |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from mouse Dystrobrevin alpha           |
| <b>Target</b>              | Dystrobrevin alpha                                                               |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| <b>Form/Appearance</b>     | Liquid                                                                           |
| <b>Concentration</b>       | 1mg/ml                                                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q9Y4J8</a>                                                           |
| <b>Dilution Range</b>      | IF=1:100-500,ICC=1:100                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |