

Product Datasheet

PPP1R12C antibody (Biotin) (orb46558)

Description	PPP1R12C antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Protein phosphatase 1 regulatory subunit 12C protein (317-450AA)
Target	PPP1R12C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZL4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.