



## Product Datasheet

### METTL14 antibody (orb46595)


|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | METTL14 antibody                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                              |
| <b>Reactivity</b>          | Human                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                        |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                      |
| <b>Immunogen</b>           | Recombinant Human N6-adenosine-methyltransferase subunit METTL14 protein (24-206AA) |
| <b>Target</b>              | METTL14                                                                             |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                         |
| <b>Form/Appearance</b>     | Liquid                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                       |
| <b>Note</b>                | For research use only                                                               |
| <b>Isotype</b>             | IgG                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                          |
| <b>Uniprot ID</b>          | <a href="#">Q9HCE5</a>                                                              |
| <b>Dilution Range</b>      | WB:1:1000-5000, IHC-P:1:20-200                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                     |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of SH-SY5Y whole ...



Immunohistochemistry of paraffin-embedde...