

Product Datasheet

HDAC1 antibody (FITC) (orb466029)

Description	HDAC1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human HDAC1
Target	HDAC1 (Nuclear Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13547
Dilution Range	Flow-Cyt=1:50-300,IF=1:100-500
Expiration Date	12 months from date of receipt.