

Product Datasheet

Tubulin beta 3 (Neuronal Marker) antibody (FITC) (orb466084)

Description	Tubulin beta 3 (Neuronal Marker) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human beta III Tubulin
Target	TUBB3 (Neuronal Marker)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13509
Dilution Range	Flow-Cyt=1µg/Test,IF=1:200-800,ICC=1:100
Expiration Date	12 months from date of receipt.