


## Product Datasheet

### ABCG1 antibody (orb46648)

|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | ABCG1 antibody                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                          |
| <b>Reactivity</b>          | Human                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ELISA, IHC                                                                      |
| <b>Immunogen</b>           | Recombinant Human ATP-binding cassette sub-family G member 1 protein (84-222AA) |
| <b>Target</b>              | ABCG1                                                                           |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                     |
| <b>Form/Appearance</b>     | Liquid                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                   |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | IgG                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                      |
| <b>Uniprot ID</b>          | <b>P45844</b>                                                                   |
| <b>Dilution Range</b>      | IHC-P:1:20-200                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



Immunohistochemistry  
of paraffin-embedde...