

Product Datasheet

beta 1 Adrenergic Receptor antibody (Biotin) (orb46663)

Description	beta 1 Adrenergic Receptor antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Beta-1 adrenergic receptor protein (401- 477AA)
Target	ADRB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08588
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.