



Product Datasheet

ALDH1A3 antibody (orb46672)

Description	ALDH1A3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Aldehyde dehydrogenase family 1 member A3 protein (2-210AA)
Target	ALDH1A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P47895
Dilution Range	WB:1:2000-10000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of U87 whole cell ...



Immunofluorescence analysis of PC-3 cell...