

## Product Datasheet

### ARHGDIB antibody (Biotin) (orb46683)

|                            |                                                                     |
|----------------------------|---------------------------------------------------------------------|
| <b>Description</b>         | ARHGDIB antibody (Biotin)                                           |
| <b>Species/Host</b>        | Rabbit                                                              |
| <b>Reactivity</b>          | Human                                                               |
| <b>Conjugation</b>         | Biotin                                                              |
| <b>Tested Applications</b> | ELISA                                                               |
| <b>Immunogen</b>           | Recombinant Human Rho GDP-dissociation inhibitor 2 protein (2-86AA) |
| <b>Target</b>              | ARHGDIB                                                             |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol      |
| <b>Form/Appearance</b>     | Liquid                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.       |
| <b>Note</b>                | For research use only                                               |
| <b>Isotype</b>             | IgG                                                                 |
| <b>Clonality</b>           | Polyclonal                                                          |
| <b>Uniprot ID</b>          | <b>P52566</b>                                                       |
| <b>Dilution Range</b>      | ELISA                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                     |